

AED Plus Trainer²

Practice an AED Plus Rescue with
Total Instructor Control



An accurate, well-controlled, and simple simulation of what it's like to use an AED Plus[®] in a real rescue.

AED Plus Trainer2 includes:

- AED Plus Trainer2 training AED
- Wireless remote controller
- One set of AED Plus Training **CPR-D•padz**[®]
- One set of AED Plus Trainer **CPR-D•padz** reusable gels
- 4 D cell batteries
- 2 AA cell batteries
- Operator's Guide

Complete Simulation Control at the Instructor's Fingertips

The ZOLL[®] AED Plus Trainer2 allows CPR and first aid instructors to have students experience what it's like to save a life with the AED Plus. Choose from four different fixed rescue scenarios or have complete manual control of the entire rescue. A wireless remote controller enables the control of multiple AED Plus Trainer2s.

With this simple remote controller, instructors can set up a rescue scenario and download to the AED Plus Trainer2 in seconds. The student pushes the on-off button on the Trainer2 unit to begin. The instructor then determines whether to shock on the next simulated heart analysis or not, whether the pads have been properly attached or not, and when the Trainer2 unit should say "PUSH HARDER" and "GOOD COMPRESSIONS."

A simple pause button lets the instructor stop the rescue altogether, confer with the student, and then pick up exactly where the simulated rescue left off.

Protocols conform to International Liaison Committee on Resuscitation (ILCOR) standards including American Heart Association (AHA) and European Resuscitation Council (ERC) 2005 Guidelines.

Trainer 2 easily switches from AHA to ERC 2005 Guidelines



AED Plus Trainer2 Specifications

Dimensions: 5.25" H x 9.5" W x 11.5" D (13.3cm H x 24.1cm W x 29.2cm D)

Weight: 4 lbs 14 oz (2.2 kg)

Power: 4 D cells (included) for Training AED (Good for more than 100 hours of training)

2 AA cells (included) for remote controller

AC adapter not available—battery operated only.

Zoll Medical
Your Safety Company
111 Quarry Road
Chambersburg PA 17202
phone 800-548-3920
fax 717-267-2369

